


SERUM

Serum derived from Atlantic salmon (*Salmo salar*)



FEATURES

- Sterile filtered (0.2 μ m), sterility tested
- Stringent collection techniques
- Antibiotic and hormone free source
- Alternative to bovine and other mammalian-derived options
- Sourced from healthy stock, tested for pathogens of regulatory concern
- Low hemoglobin content

APPLICATION

- Culture of fish cells, compatible with insect cell lines
- Cell culture research

SPECIFICATIONS

- Frozen liquid
- 0.2 μ m filtered
- Available in 50 ml and 100 ml aliquots
- Endotoxin & hemoglobin screening available
- Certificate of Analysis provided
- Custom preparations to meet different requirements

PRODUCT NO.	015	STORAGE	$\leq -60^{\circ}\text{C}$
--------------------	-----	----------------	----------------------------

UNIT OF SALE	50 ml bottle	SHELF LIFE	2 years
---------------------	--------------	-------------------	---------