


SALMON PLASMA

Plasma derived from Atlantic salmon (*Salmo salar*)



FEATURES

- Sterile filtered (0.2µm), sterility tested
- Antibiotic and hormone free source
- Alternative to bovine and other mammalian-derived options
- Sourced from healthy stock, tested for pathogens of regulatory concern
- Available with citrate, heparin, and EDTA

APPLICATION For research use only

FORMULATION Frozen liquid

- 0.2 µm filtered
- Available in 50 ml aliquots; Larger volumes and custom formulations are available upon request
- Variety of anticoagulants:
 - (028) sodium citrate
 - (030) sodium EDTA
 - (032) dipotassium EDTA
 - (034) tripotassium EDTA
 - (036) sodium heparin

PRODUCT NO. 028, 030, 032, 034, 036

STORAGE ≤-60°C

UNIT OF SALE 50 ml bottle

SHELF LIFE 2 years from manufacture